

Labycare-Elbow Crutches (SKSS-RP480)



Elbow Crutches Manufacturer by LABYCARE

Description :

Elbow crutches are manufactured from aluminium.

- * It has a cuff at the top to go around the forearm.
- * Elbow crutches are used by slipping your arm into the cuff and holding the hand grip.
- * The cuff is made of plastic and can be closed or a half-circle with a opening in the front, allowing the forearm to slip out in case of a fall.

We are Manufacturer, Exporters and Suppliers of Elbow Crutches at wholesale competitive prices from India.

Specifications :

Elbow crutches are manufactured from aluminium.

- * It has a cuff at the top to go around the forearm.
- * Elbow crutches are used by slipping your arm into the cuff and holding the hand grip.
- * The cuff is made of plastic and can be closed or a half-circle with a opening in the front, allowing the forearm to slip out in case of a fall.

We are Manufacturer, Exporters and Suppliers of Elbow Crutches at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com