

Labycare-Histamine Chamber (SKSSHC913)



Histamine Chamber

Description :

Histamine Chamber

Specifications : Histamine Chamber SKSS Histamine Chamber for Guinea Pigs Large chamber with special perforated removable partition, with inbuilt pressureised filtered air delivery electric pump and fine mist nebulizer eliminates need for oxygen cylinder manual efforts. Size 14x12x9 inch.

TECHNICAL SPECIFICATION:

Shape	Rectang
Automation Type	Manual
Brand	SKSS
Material	Glass
Accuracy (millimetre)	+ -2%
Usage	Laborato

labycare

www.labycare.com | Email: labycare@gmail.com