


Labycare- Mammalian heart assembly (SK-DA)



Digital Assembly for perfusion of excised Mammalian heart / Langendroff's assembly - Single unit

Description :

Specifications :	Digital Assembly forSKSS Digital Assembly for perfusion of excised perfusion of excisedMammalian heart with stainless steel fittings, Mammalian heart /perspex bath with PID Digital control,(+/-0.5*C) Langendroff's assemblyelectrical heating unit with double reservoir - Single unit unit. SINGLE UNIT.
Brand	SKSS

labycare

www.labycare.com | Email: labycare@gmail.com