

## Labycare-Frog Splint (SKSS-RP537)



Frog Splint Manufacturer by LABYCARE

### **Description :**

Frog Splint is recommended for hyper-flexion injuries & collateral ligament injuries of Distal Inter Phalangeal joint.

Frog Splint supports and stabilizes the finger.

Frog Splint is made of malleable aluminium with upward curve at tip, protects the distal end of finger. Adjustable velcro strap. Easy application.

Foam padding ensures comfortable fit.

Common for left & right hand finger.

We are Exporters and Suppliers of Frog Splint at wholesale competitive prices from India.

### **Specifications :**

Frog Splint is recommended for hyper-flexion injuries & collateral ligament injuries of Distal Inter Phalangeal joint.

Frog Splint supports and stabilizes the finger.

Frog Splint is made of malleable aluminium with upward curve at tip, protects the distal end of finger. Adjustable velcro strap. Easy application.

Foam padding ensures comfortable fit.

Common for left & right hand finger.

We are Exporters and Suppliers of Frog Splint at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)