

## Labycare-Extended Baseball Splint (SKSS-RP540)



Extended Baseball Splint Manufacturer by LABYCARE

### **Description :**

Extended Baseball Splint is recommended for finger fracture, ligament injury, finger deformity .

Supports, protects and immobilizes the finger and metacarpophalangeal joint.

Maintains interphalangeal joint in correct position. Made of malleable aluminium splint which can be shaped to maintain MCP joint in flexed position.

Easy to wear with comfortable foam padding.

Common for left & right hand finger.

We are Exporters and Suppliers of Extended Baseball Splint at wholesale competitive prices from India.

### **Specifications :**

Extended Baseball Splint is recommended for finger fracture, ligament injury, finger deformity .

Supports, protects and immobilizes the finger and metacarpophalangeal joint.

Maintains interphalangeal joint in correct position. Made of malleable aluminium splint which can be shaped to maintain MCP joint in flexed position.

Easy to wear with comfortable foam padding.

Common for left & right hand finger.

We are Exporters and Suppliers of Extended Baseball Splint at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)