

Labycare-Stax Mallet (SKSS-RP541)



Stax Mallet Manufacturer by LABYCARE

Description :

Stax Mallet is useful in finger deformity.

Stax Mallet supports the distal interphalangeal joint in hyper extension, keeping proximal interphalangeal joint free.

Stax Mallet is made of engineering plastic. Well ventilated holes keeps the skin cool and dry. Easy application and removal.

It keeps the finger in 10 degree of hyper extension. Pre- cut design for perfect fit.

We are Exporters and Suppliers of Stax Mallet at wholesale competitive prices from India.

Specifications :

Stax Mallet is useful in finger deformity.

Stax Mallet supports the distal interphalangeal joint in hyper extension, keeping proximal interphalangeal joint free.

Stax Mallet is made of engineering plastic. Well ventilated holes keeps the skin cool and dry. Easy application and removal.

It keeps the finger in 10 degree of hyper extension. Pre- cut design for perfect fit.

We are Exporters and Suppliers of Stax Mallet at wholesale competitive prices from India.

