

Labycare-long arm splint (SKSS-RP543)



long arm splint Manufacturer by labycare

Description :

Long Arm Splint is used immediately post-injury or in emergency splinting of fracture / dislocation of upperlimb.

Long Arm Splint supports, protects and immobilizes affected area.

Long Arm Splint is made of malleable metal to shape the splint according to need.

We are Exporters and Suppliers of Long Arm Splint at wholesale competitive prices from India.

Specifications :

Long Arm Splint is used immediately post-injury or in emergency splinting of fracture / dislocation of upperlimb.

Long Arm Splint supports, protects and immobilizes affected area.

Long Arm Splint is made of malleable metal to shape the splint according to need.

We are Exporters and Suppliers of Long Arm Splint at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com