

Labycare-Foot Drop Splint (SKSS-RP547)



Foot Drop Splint Manufacturer by LABYCARE

Description :

Foot Drop Splint is used for foot drop condition, ankle or plantar flexion contracture, post-operative care, burns or in peripheral nerve damage.

Foot Drop Splint supports the foot during walking by preventing foot slap at heel strike and assist toe clearance in swing phase.

Semi rigid construction with hook and loop closures.

Should be used with footwear. Easy application.

Helps in walking. Can be worn with or without shoe. Easy to use.

Inbuilt hook loop closures & padding ensure better fitting.

Different for left & right Leg.

We are Exporters and Suppliers of Foot Drop Splint at wholesale competitive prices from India.

Specifications :

Foot Drop Splint is used for foot drop condition, ankle or plantar flexion contracture, post-operative care, burns or in peripheral nerve damage.

Foot Drop Splint supports the foot during walking by preventing foot slap at heel strike and assist toe clearance in swing phase.

Semi rigid construction with hook and loop closures.

Should be used with footwear. Easy application.

Helps in walking. Can be worn with or without shoe. Easy to use.

Inbuilt hook loop closures & padding ensure better fitting.

Different for left & right Leg.

We are Exporters and Suppliers of Foot Drop Splint at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com