

Labycare-Cervical Traction Kit for Sitting (SKSS-RP548)



Cervical Traction Kit for Sitting Manufacturer by LABYCARE

Description :

Cervical Traction Kit is recommended for sprain or dislocation in cervical region.

It provides traction to the cervical & upper thoracic vertebrae.

Cervical Traction Kit relieves pressure by stretching musculature & connective tissue of cervical region.

Cervical Traction Kit consists of soft padded head holder, spreader bar, pulley set, 3 meters long cord, water weight bag & door stand.

We are Exporters and Suppliers of Cervical Traction Kit for Sitting at wholesale competitive prices from India.

Specifications :

Cervical Traction Kit is recommended for sprain or dislocation in cervical region.

It provides traction to the cervical & upper thoracic vertebrae.

Cervical Traction Kit relieves pressure by stretching musculature & connective tissue of cervical region.

Cervical Traction Kit consists of soft padded head holder, spreader bar, pulley set, 3 meters long cord, water weight bag & door stand.

We are Exporters and Suppliers of Cervical Traction Kit for Sitting at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com