

Labycare-Skin & Foot Traction Kit (SKSS-RP561)



Skin & Foot Traction Kit Manufacturer by LABYCARE

Description :

Skin & Foot Traction Kit provides traction.

Stripes provide maximum aeration to the skin. Elastic cloth is fastened with low pressure.

Skin & Foot Traction Kit consists of water weight bag, 3 meters long cord, traction pulley bracket & skin foot belt.

We are Exporters and Suppliers of Skin & Foot Traction Kit at wholesale competitive prices from India.

Specifications :

Skin & Foot Traction Kit provides traction.

Stripes provide maximum aeration to the skin. Elastic cloth is fastened with low pressure.

Skin & Foot Traction Kit consists of water weight bag, 3 meters long cord, traction pulley bracket & skin foot belt.

We are Exporters and Suppliers of Skin & Foot Traction Kit at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com