

Labycare-adhesive skin traction kit (SKSS-RP551)



adhesive skin traction kit manufacturer by labycare

Description :

Adhesive Skin Traction Kit is most effective method of skin traction. Pressure sensitive acrylic adhesive belt to give maximum adhesion with the minimum allergy related complications.

Adhesive Skin Traction Kit consists of water weight bag. 3 meter long cord, traction pulley bracket & adhesive skin foot belt.

We are Exporters and Suppliers of Adhesive Skin Traction Kit at wholesale competitive prices from India.

Specifications :

Adhesive Skin Traction Kit is most effective method of skin traction. Pressure sensitive acrylic adhesive belt to give maximum adhesion with the minimum allergy related complications.

Adhesive Skin Traction Kit consists of water weight bag. 3 meter long cord, traction pulley bracket & adhesive skin foot belt.

We are Exporters and Suppliers of Adhesive Skin Traction Kit at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com