

Labycare-leg traction kit (SKSS-RP552)



leg traction kit manufacturer by labycare

Description :

Leg Traction Kit provides comfort with good skin traction.

Leg traction brace has advantage of easy application & removal.

It is applied below knee traction with very large holding surface area.

Leg Traction Kit consists of leg traction brace, 3 meters long cord, traction pulley bracket & water weight bag.

We are Exporters and Suppliers of Leg Traction Kit at wholesale competitive prices from India.

Specifications :

Leg Traction Kit provides comfort with good skin traction.

Leg traction brace has advantage of easy application & removal.

It is applied below knee traction with very large holding surface area.

Leg Traction Kit consists of leg traction brace, 3 meters long cord, traction pulley bracket & water weight bag.

We are Exporters and Suppliers of Leg Traction Kit at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com