

Labycare-Pelvic Traction (SKSS-PR555)



Pelvic Traction Manufacturer by LABYCARE

Description :

Pelvic Traction includes pelvic belt, spreader bar & 3 meters long cord.

Pelvic Traction provides traction to lumbo - sacral & lower thoracic regions.

Pelvic Traction provides uniform & controlled traction.

Pelvic Traction is made of coated fabric with foam padding, lining of sweat absorbing fabric for comfort. Velcro closure for easy application.

We are Exporters and Suppliers of Pelvic Traction at wholesale competitive prices from India.

Specifications :

Pelvic Traction includes pelvic belt, spreader bar & 3 meters long cord.

Pelvic Traction provides traction to lumbo - sacral & lower thoracic regions.

Pelvic Traction provides uniform & controlled traction.

Pelvic Traction is made of coated fabric with foam padding, lining of sweat absorbing fabric for comfort. Velcro closure for easy application.

We are Exporters and Suppliers of Pelvic Traction at wholesale competitive prices from India.