

Labycare-Tripod Walking Sticks (SKSS-RP567)



Tripod Walking Sticks Manufacturer by LABYCARE

Description :

Tripod Walking Sticks made of Aluminium.

- * Tripod base provides greater stability.
- * Angled neck reduces strain on wrists.
- * Non-marking, slip resistant rubber tips.
- * Height adjustment 79 cm to 96 cm.
- * Base measurement 13 x 13 x 13 cm..

We are Manufacturer, Exporters and Suppliers of Tripod Walking Sticks at wholesale competitive prices from India.

Specifications :

Tripod Walking Sticks made of Aluminium.

- * Tripod base provides greater stability.
- * Angled neck reduces strain on wrists.
- * Non-marking, slip resistant rubber tips.
- * Height adjustment 79 cm to 96 cm.
- * Base measurement 13 x 13 x 13 cm..

We are Manufacturer, Exporters and Suppliers of Tripod Walking Sticks at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com