

Labycare-Quadripod Walking Sticks Standard (SKSS-RP568)



Quadripod Walking Sticks Standard

Description :

Quadripod Walking Sticks Standard is made of Aluminium.

- * Quadruped base provides greater stability.
- * Swan neck reduces strain on wrists.
- * Casting aluminium base.
- * Adjustable height.
- * Non-marking, slip resistant rubber tips.
- * Max User weight 125 kgs.
- * Height adjustment 72cm to 95cm.

We are Manufacturer, Exporters and Suppliers of Quadripod Walking Sticks Standard at wholesale competitive prices from India.

Specifications :

Quadripod Walking Sticks Standard is made of Aluminium.

- * Quadruped base provides greater stability.
- * Swan neck reduces strain on wrists.
- * Casting aluminium base.

- * Adjustable height.
- * Non-marking, slip resistant rubber tips.
- * Max User weight 125 kgs.
- * Height adjustment 72cm to 95cm.

We are Manufacturer, Exporters and Suppliers of Quadripod Walking Sticks Standard at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com