

Labycare-Wheelchair Standard Powder Coated (SKSS-RP571)



Wheelchair Standard Powder Coated Manufacturer by LABYCARE

Description :

Wheelchair Standard PC is made of Mild Steel Frame

Mild Steel Frame is further Chemical Treated Epoxy Powder Coated

Fixed PVC Armrest is fitted with Wheelchair

PVC Seat and Back is fitted.

Fixed Aluminium Footrest is fitted to the wheelchair

Regular Brakes are available on the rear wheels.

We are Manufacturer, Exporters and Suppliers of Wheelchair Standard PC at wholesale competitive prices from India

Specifications :

Wheelchair Standard PC is made of Mild Steel Frame

Mild Steel Frame is further Chemical Treated Epoxy Powder Coated

Fixed PVC Armrest is fitted with Wheelchair

PVC Seat and Back is fitted.

Fixed Aluminium Footrest is fitted to the wheelchair

Regular Brakes are available on the rear wheels.

We are Manufacturer, Exporters and Suppliers of Wheelchair Standard PC at wholesale competitive prices from India

labycare

www.labycare.com | Email: labycare@gmail.com