

Labycare-Wheelchair with Handbrake (SKSS982A)



Wheelchair with Handbrake Manufacturer by **LABYCARE** in India.

Description :

Labycare Wheelchair _____

Specifications :

- Wheelchair Manual with Hand Brakes is made of Thick Mild Steel Frame.
- Mild Steel Frame is further Chrome Plated.
- Fixed PVC Armrest is fitted with Wheelchair.
- Soft Cushioned PU Seat & Back with pocket is fitted.
- Plastic Hand Rim is fitted on both the rear wheels.
- Adjustable Aluminium Footrest is fitted to the wheelchair.
- Brakes are available on Hand Grip for the attendant and rear wheels for the patient.

1. Front Castor Size: 20 cm
2. Rear Wheel Size: 60 cm
3. Seat Depth: 44 cm
4. Seat Width: 50 cm
5. Seat Height: 46 cm
6. Loading Weight: 100 kg

