

Labycare-Backrest Reclining Wheelchair (SKSS-RP577)



Backrest Reclining Wheelchair Manufacturer by LABYCARE

Description :

Backrest Reclining Wheelchair is made of Mild Steel Frame

Wheelchair Backrest Reclining 90 -180 degree

Mild Steel Frame is further Chrome Plated

Fixed PVC Armrest is fitted with Wheelchair

PVC Seat and Back is fitted

Footrest are Detachable & Elevating

Aluminium Footrest is fitted to the wheelchair

United Brakes are available

We are Manufacturer, Exporters and Suppliers of Backrest Reclining Wheelchair at wholesale competitive prices from India

Specifications :

Backrest Reclining Wheelchair is made of Mild Steel Frame

Wheelchair Backrest Reclining 90 -180 degree

Mild Steel Frame is further Chrome Plated

Fixed PVC Armrest is fitted with Wheelchair

PVC Seat and Back is fitted

Footrest are Detachable & Elevating

Aluminium Footrest is fitted to the wheelchair

United Brakes are available

We are Manufacturer, Exporters and Suppliers of Backrest Reclining Wheelchair at wholesale competitive prices from India

labycare

www.labycare.com | Email: labycare@gmail.com