

Labycare-Instrument Cleaner (SKSSIC226)



Instrument Cleaner

Description : • Base on H₂O, peroxide, enzymes, and chelating agents. It is suitable for cleaning surgical instruments, glassware, plastic, metal, and rubber goods.

Specifications :

- Combination of enzymes for an enhanced cleaning action
- Drastically removes organic residues
- Very economical: Up to 0,5% dilution
- Suitable for the cleaning in ultrasonic bath, in instrument washers, in immersion bath, by conventional hand washing as well as in automated endoscope disinfectors
- Does not foam, prevents corrosion and instrument discoloration
- Appropriate for heat-resistant as well as heat-sensitive instruments
- 8oz/2 liter bottles. neutral pH formula for use in ultrasonic or manual cleaning.
- Designed for use on surgical instruments, glassware, plastic, metal or rubber goods.
- One bottle yields 32 gallons/121 liters of solution.
- Case of 12.

labycare

www.labycare.com | Email: labycare@gmail.com