

Labycare-Aluminium Wheelchair Heavy Duty (SKSS-RP578)



Aluminium Wheelchair Heavy Duty Manufacturer by labycare

Description :

Wheelchair is made of Heavy Duty Aluminium Frame.

Aluminium Frame is further treated to Power Coating Frame Colour

Flipup PU Armrest is fitted with Wheelchair

Nylon Seat and Back is fitted.

Detachable Footrest is fitted to the wheelchair

Regular Brakes are available.

We are Manufacturer, Exporters and Suppliers of Aluminium Wheelchair Heavy Duty at wholesale competitive prices from India

Specifications :

Wheelchair is made of Heavy Duty Aluminium Frame.

Aluminium Frame is further treated to Power Coating Frame Colour

Flipup PU Armrest is fitted with Wheelchair

Nylon Seat and Back is fitted.

Detachable Footrest is fitted to the wheelchair

Regular Brakes are available.

We are Manufacturer, Exporters and Suppliers of Aluminium Wheelchair Heavy Duty at wholesale competitive prices from India

labycare

www.labycare.com | Email: labycare@gmail.com