

Labycare-Wheelchair Electric (SKSS-RP581)



Wheelchair Electric Manufacturer by LABYCARE

Description :

Electric Wheel Chair is made of Mild Steel Frame

Mild Steel Frame is further Chemical Treated Epoxy Powder Coated

Flip up Armrest is fitted with Wheelchair

Detachable Footrest is fitted to the wheelchair

Drop Back Handle is available

Castors 200 x 50 PU

Drivewheel 12 inch PU

We are Manufacturer, Exporters and Suppliers of Wheelchair Electric at wholesale competitive prices from India

Specifications :

Electric Wheel Chair is made of Mild Steel Frame

Mild Steel Frame is further Chemical Treated Epoxy Powder Coated

Flip up Armrest is fitted with Wheelchair

Detachable Footrest is fitted to the wheelchair

Drop Back Handle is available

Castors 200 x 50 PU

Drivewheel 12 inch PU

We are Manufacturer, Exporters and Suppliers of Wheelchair Electric at wholesale competitive prices from India

labycare

www.labycare.com | Email: labycare@gmail.com