

## Labycare-High Vacuum Suction Machine (SKSS-SM03)



High Vacuum Suction Machine Manufacturer by LABYCARE

### **Description :**

Suction Machine StanVac is High Vacuum machine.

Suction Machine Cabinet is made by mild steel with electrostatically powder coated Anti corrosive mounted on 63 mm diameter on (4) castors

Suction Machine StanVac is fitted with Reuseable Filter Jars : 2 x 2.5 Ltrs.

The Polycarbonate Jars are fitted with lids with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 100 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Twin Headed Piston Pump

Vacuum Capacity is -730 mm Hg  $\pm$  10% at 95-100 LPM which is controlled by knob.

Noise Level < 50 dB A  $\pm$  3 Almost Whispers

Motor is thermally protected

Power : 220/230V AC , 50 Hz , 480W

The valve is made of Stainless Steel.

We are Manufacturer, Exporters and Suppliers of High Vacuum Suction Machine StanVac at wholesale competitive prices from India.

**Specifications :**

Suction Machine StanVac is High Vacuum machine.

Suction Machine Cabinet is made by mild steel with electrostatically powder coated Anti corrosive mounted on 63 mm diameter on (4) castors

Suction Machine StanVac is fitted with Reuseable Filter Jars : 2 x 2.5 Ltrs.

The Polycarbonate Jars are fitted with lids with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 100 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Twin Headed Piston Pump

Vacuum Capacity is -730 mm Hg  $\pm$  10% at 95-100 LPM which is controlled by knob.

Noise Level < 50 dB A  $\pm$  3 Almost Whispers

Motor is thermally protected

Power : 220/230V AC , 50 Hz , 480W

The valve is made of Stainless Steel.

We are Manufacturer, Exporters and Suppliers of High Vacuum Suction Machine StanVac at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)