

Labycare-High Vacuum Suction Machine DeVac 2 (SKSS-SM06)



High Vacuum Suction Machine DeVac 2 Manufacturer by labycare

Description :

Suction Machine DeVac 2 is High Vacuum machine.

Suction Machine Cabinet is made by mild steel with electrostatically powder coated Anti corrosive mounted on 63 mm diameter on (4) castors

Suction Machine DeVac 2 is fitted with Reuseable Filter Jars : 2 x 4 Ltrs.

The Polycarbonate Jars are fitted with lids with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 63 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Twin Headed Piston Pump

Vacuum Capacity is -740 mm Hg \pm 10% at 120-125 LPM which is controlled by knob.

Noise Level < 50 dB A \pm 3 Almost Whispers

Motor is Thermally protected

Power : 220/230V AC , 50 Hz , 820W

Valves are Stainless Steel

We are Manufacturer, Exporters and Suppliers of High Vacuum Suction Machine DeVac 2 at wholesale competitive prices from India.

Specifications :

Suction Machine DeVac 2 is High Vacuum machine.

Suction Machine Cabinet is made by mild steel with electrostatically powder coated Anti corrosive mounted on 63 mm diameter on (4) castors

Suction Machine DeVac 2 is fitted with Reuseable Filter Jars : 2 x 4 Ltrs.

The Polycarbonate Jars are fitted with lids with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 63 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Twin Headed Piston Pump

Vacuum Capacity is -740 mm Hg \pm 10% at 120-125 LPM which is controlled by knob.

Noise Level < 50 dB A \pm 3 Almost Whispers

Motor is Thermally protected

Power : 220/230V AC , 50 Hz , 820W

Valves are Stainless Steel

We are Manufacturer, Exporters and Suppliers of High Vacuum Suction Machine DeVac 2 at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com