

Labycare-Electric Suction Machine MS StanVac (SKSS-SM08)



Electric Suction Machine MS StanVac Manufacturer by LABYCARE

Description :

Suction Machine Cabinet is made by mild steel with electrostatically powder coated mounted on 50 mm diameter on (4) Rubber castors

Suction Machine StanVac has stainless steel top.

Suction Machine StanVac is fitted with Reuseable Filter Jars : 2 x 1.5 Ltrs.

The Polycarbonate (PC) Jars are fitted with synthetic rubber lids with over flow safety device.

This Suction machine is Oil Immersed Motorised Noiseless Vacuum Pump

Vacuum Gauge has 63 mm diameter graduated in mm Hg.

Vacuum is -700 mm Hg \pm 10% which is controlled by knob.

Power : 230V-50 Hz

Foot Switch Socket is provided for foot switch.

We are Manufacturer, Exporters and Suppliers of Electric Suction Machine MS StanVac at wholesale competitive prices from India.

Specifications :

Suction Machine Cabinet is made by mild steel with electrostatically powder coated mounted on 50 mm diameter on (4) Rubber castors

Suction Machine StanVac has stainless steel top.

Suction Machine StanVac is fitted with Reuseable Filter Jars : 2 x 1.5 Ltrs.

The Polycarbonate (PC) Jars are fitted with synthetic rubber lids with over flow safety device.

This Suction machine is Oil Immersed Motorised Noiseless Vacuum Pump

Vacuum Gauge has 63 mm diameter graduated in mm Hg.

Vacuum is -700 mm Hg \pm 10% which is controlled by knob.

Power : 230V-50 Hz

Foot Switch Socket is provided for foot switch.

We are Manufacturer, Exporters and Suppliers of Electric Suction Machine MS StanVac at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com