

## Labycare-Electric Suction Machine ABS Body (SKSS-SM14)



Electric Suction Machine ABS Body Manufacturer by LABYCARE

### **Description :**

Suction Machine Cabinet is made by Acrylonitrile-Butadiene-Styrene (ABS) Plastic

Anti corrosive mounted on 63 mm diameter on (4) castors

Suction Machine ABS Body is fitted with Reuseable Filter Jars : 2 x 2 Ltrs.

The Polycarbonate (PC) Jars are fitted with lids with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 75 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Twin Headed Piston Pump

Vacuum Capacity is -710 mm Hg  $\pm$  10% at 32-35 LPM which is controlled by knob.

Noise Level < 50 dB A  $\pm$  3 Almost Whispers

Motor is thermally protected

Power : 220/230V AC , 50 Hz , 130W

The Valve is Stainless Steel

We are Manufacturer, Exporters and Suppliers of Electric Suction Machine ABS Body at wholesale competitive prices from India.

**Specifications :**

Suction Machine Cabinet is made by Acrylonitrile-Butadiene-Styrene (ABS) Plastic

Anti corrosive mounted on 63 mm diameter on (4) castors

Suction Machine ABS Body is fitted with Reuseable Filter Jars : 2 x 2 Ltrs.

The Polycarbonate (PC) Jars are fitted with lids with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 75 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Twin Headed Piston Pump

Vacuum Capacity is -710 mm Hg  $\pm$  10% at 32-35 LPM which is controlled by knob.

Noise Level < 50 dB A  $\pm$  3 Almost Whispers

Motor is thermally protected

Power : 220/230V AC , 50 Hz , 130W

The Valve is Stainless Steel

We are Manufacturer, Exporters and Suppliers of Electric Suction Machine ABS Body at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)