

## Labycare-Electric Suction Machine Stainless Steel (SKSS-SM15)



Electric Suction Machine Stainless Steel Manufacturer by LABYCARE

### **Description :**

Suction Machine complete housing is made by Stainless Steel

Anti corrosive mounted on 63 mm diameter on (4) castors

Suction Machine ABS Body is fitted with Reuseable Filter Jars : 2 x 1.5 Ltrs.

The Polycarbonate (PC) Jars are fitted with lids with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 75 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Immersed Rotary Vane

Vacuum Capacity is -710 mm Hg  $\pm$  10% at 32-35 LPM which is controlled by knob.

Noise Level < 50 dB A  $\pm$  3 Almost Whispers

Motor is 1/4 HP FHP

Power : 220/230V AC , 50 Hz , 180W

The RPM is 1440

We are Manufacturer, Exporters and Suppliers of Electric Suction Machine Stainless Steel at wholesale competitive prices from India.

**Specifications :**

Suction Machine complete housing is made by Stainless Steel

Anti corrosive mounted on 63 mm diameter on (4) castors

Suction Machine ABS Body is fitted with Reuseable Filter Jars : 2 x 1.5 Ltrs.

The Polycarbonate (PC) Jars are fitted with lids with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 75 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Immersed Rotary Vane

Vacuum Capacity is -710 mm Hg  $\pm$  10% at 32-35 LPM which is controlled by knob.

Noise Level < 50 dB A  $\pm$  3 Almost Whispers

Motor is 1/4 HP FHP

Power : 220/230V AC , 50 Hz , 180W

The RPM is 1440

We are Manufacturer, Exporters and Suppliers of Electric Suction Machine Stainless Steel at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)