

Labycare-Suction Machine Stand Model Stanvac 1 (SKSS-SM16)



Suction Machine Stand Model Stanvac 1 Manufacturer by LABYCARE

Description :

Suction Machine Stand Model Stanvac 1 is made of by mild steel with electrostatically powder coated Steel

The Unit is mounted on Mild Steel Stand with 50 mm diameter on (4) Rubber castors

Suction Machine Stanvac 1 has Oil Rotary Noiseless Pump

Suction Machine Stanvac 1 is fitted with Reuseable Filter Jars : 2 x 2 Ltrs.

The Polycarbonate Jars are fitted on the top of the Mild Steel Stand base with lids with over flow safety device.

This machine is versatile for Surgical Applications

Vacuum Gauge has 63 mm diameter graduated in mm Hg.

Vacuum is -700 mm Hg \pm 10% which is controlled by knob.

Power : 230V-50 Hz

We are Manufacturer, Exporters and Suppliers of Suction Machine Stand Model Stanvac 1 at wholesale competitive prices from India.

Specifications :

Suction Machine Stand Model Stanvac 1 is made of by mild steel with electrostatically powder coated

Steel

The Unit is mounted on Mild Steel Stand with 50 mm diameter on (4) Rubber castors

Suction Machine Stanvac 1 has Oil Rotary Noiseless Pump

Suction Machine Stanvac 1 is fitted with Reuseable Filter Jars : 2 x 2 Ltrs.

The Polycarbonate Jars are fitted on the top of the Mild Steel Stand base with lids with over flow safety device.

This machine is versatile for Surgical Applications

Vacuum Gauge has 63 mm diameter graduated in mm Hg.

Vacuum is -700 mm Hg \pm 10% which is controlled by knob.

Power : 230V-50 Hz

We are Manufacturer, Exporters and Suppliers of Suction Machine Stand Model Stanvac 1 at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com