

Labycare-Suction Machine Stand Model EcoVac 1 (SKSS-SM17)



Suction Machine Stand Model EcoVac 1 Manufacturer by LABYCARE

Description :

Suction Machine Stand Model EcoVac 1 is made of ABS Plastic mounted on Mild Steel Stand with 50 mm diameter on (5) Rubber castors

Suction Machine EcoVac 1 has Oil Free Diaphragm Pump not to be lubricated

Suction Machine EcoVac 1 is fitted with Reuseable Filter Jars : 2 Jars.

The Polycarbonate Jars are fitted on the sides of the Metal Steel stand with lids with over flow safety device.

Vacuum Gauge has 50 mm diameter graduated in mm Hg.

Vacuum is -600 mm Hg \pm 10% which is controlled by knob.

The Suction Capacity is 15 Ltr/min

Power : 230V-50 Hz

We are Manufacturer, Exporters and Suppliers of Suction Machine Stand Model EcoVac 1 at wholesale competitive prices from India.

Specifications :

Suction Machine Stand Model EcoVac 1 is made of ABS Plastic mounted on Mild Steel Stand with 50 mm diameter on (5) Rubber castors

Suction Machine EcoVac 1 has Oil Free Diaphragm Pump not to be lubricated

Suction Machine EcoVac 1 is fitted with Reuseable Filter Jars : 2 Jars.

The Polycarbonate Jars are fitted on the sides of the Metal Steel stand with lids with over flow safety device.

Vacuum Gauge has 50 mm diameter graduated in mm Hg.

Vacuum is -600 mm Hg \pm 10% which is controlled by knob.

The Suction Capacity is 15 Ltr/min

Power : 230V-50 Hz

We are Manufacturer, Exporters and Suppliers of Suction Machine Stand Model EcoVac 1 at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com