

## Labycare-AC/DC Portable Suction Machine EcoVac 4 (SKSS-SM21)



Portable Suction Machine EcoVac is made by LABYCARE

### **Description :**

Portable Suction Machine EcoVac is made by ABS Plastic with Stainless Steel back

Portable Suction Machine EcoVac is fitted with Reuseable Filter Jars : 1 x 500 ml.

The Polycarbonate Jar is fitted with lids on the side of the Suction Unit with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 50 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Rocker Piston Pump

Vacuum Capacity is -650 mm Hg  $\pm$  10% at 25-30 LPM which is controlled by knob.

Noise Level < 50 dB A  $\pm$  3 Almost Whispers

Motor PMDC

Power : 220/230V AC/ 12V DC

Portable Suction Machine EcoVac has a battery backup of 4 hours

We are Manufacturer, Exporters and Suppliers of AC/DC Portable Suction Machine EcoVac 4 at wholesale competitive prices from India.

### **Specifications :**

Portable Suction Machine EcoVac is made by ABS Plastic with Stainless Steel back

Portable Suction Machine EcoVac is fitted with Reuseable Filter Jars : 1 x 500 ml.

The Polycarbonate Jar is fitted with lids on the side of the Suction Unit with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 50 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Rocker Piston Pump

Vacuum Capacity is -650 mm Hg  $\pm$  10% at 25-30 LPM which is controlled by knob.

Noise Level < 50 dB A  $\pm$  3 Almost Whispers

Motor PMDC

Power : 220/230V AC/ 12V DC

Portable Suction Machine EcoVac has a battery backup of 4 hours

We are Manufacturer, Exporters and Suppliers of AC/DC Portable Suction Machine EcoVac 4 at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)