

Labycare-Portable Suction Machine StanVac (SKSS-SM22)



Portable Suction Machine StanVac Manufacturer by LABYCARE

Description :

Portable Suction Machine StanVac is made completely by ABS Plastic

Portable Suction Machine StanVac is fitted with Reuseable Filter Jars : 1 x 1.5 Ltrs.

The Polycarbonate Jar is fitted with lids on the side of the Suction Unit with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 50 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Rocker Piston Pump

Vacuum Capacity is -600 mm Hg \pm 10% at 15-20 LPM which is controlled by knob.

Noise Level < 50 dB A \pm 3 Almost Whispers

Motor Shared Pole

Power : 220/230V AC, 50HZ , 90W

The Valve is Synthetic Rubber

We are Manufacturer, Exporters and Suppliers of Portable Suction Machine StanVac at wholesale competitive prices from India.

Specifications :

Portable Suction Machine StanVac is made completely by ABS Plastic

Portable Suction Machine StanVac is fitted with Reuseable Filter Jars : 1 x 1.5 Ltrs.

The Polycarbonate Jar is fitted with lids on the side of the Suction Unit with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 50 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Rocker Piston Pump

Vacuum Capacity is -600 mm Hg \pm 10% at 15-20 LPM which is controlled by knob.

Noise Level < 50 dB A \pm 3 Almost Whispers

Motor Shared Pole

Power : 220/230V AC, 50HZ , 90W

The Valve is Synthetic Rubber

We are Manufacturer, Exporters and Suppliers of Portable Suction Machine StanVac at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com