

Labycare-Portable Slow Suction Machine EcoVac (SKSS-SM24)



Portable Slow Suction Machine EcoVac Manufacturer by LABYCARE

Description :

Portable Slow Suction Machine EcoVac is made of ABS Plastic

Portable Slow Suction Machine EcoVac has Oil Free Diaphragm Pump not to be lubricated

Portable Slow Suction Machine EcoVac is fitted with Reuseable Filter Jars : 2 Jars.

The Polycarbonate Jars are fitted on the top of the ABS Plastic cabinet with lids with over flow safety device.

Vacuum Gauge has 50 mm diameter graduated in mm Hg.

Vacuum is -600 mm Hg \pm 10% which is controlled by knob.

The Suction Capacity is 15 Ltr/min

Power : 230V-50 Hz

We are Manufacturer, Exporters and Suppliers of Portable Slow Suction Machine EcoVac at wholesale competitive prices from India.

Specifications :

Portable Slow Suction Machine EcoVac is made of ABS Plastic

Portable Slow Suction Machine EcoVac has Oil Free Diaphragm Pump not to be lubricated

Portable Slow Suction Machine EcoVac is fitted with Reuseable Filter Jars : 2 Jars.

The Polycarbonate Jars are fitted on the top of the ABS Plastic cabinet with lids with over flow safety

device.

Vacuum Gauge has 50 mm diameter graduated in mm Hg.

Vacuum is -600 mm Hg \pm 10% which is controlled by knob.

The Suction Capacity is 15 Ltr/min

Power : 230V-50 Hz

We are Manufacturer, Exporters and Suppliers of Portable Slow Suction Machine EcoVac at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com