

## Labycare-Portable Slow Suction Machine StanVac 1 (SKSS-SM25)



Portable Slow Suction Machine StanVac 1 Manufacturer by LABYCARE

### **Description :**

Portable Slow Suction Machine StanVac 1 is made of mild steel with electrostatically powder coated mounted on Rubber pad

Portable Slow Suction Machine StanVac 1 has Oil Immersed Motorised Noiseless Pump

Portable Slow Suction Machine StanVac 1 is fitted with Reuseable Filter Jars : 2 Jars.

The Polycarbonate Jars are fitted on the top of the Mild Steel Base with lids with over flow safety device.

Vacuum Gauge has 50 mm diameter graduated in mm Hg.

Vacuum is -700 mm Hg  $\pm$  10% which is controlled by knob.

The Suction Capacity is 20 Ltr/min

Power : 230V-50 Hz

We are Manufacturer, Exporters and Suppliers of Portable Slow Suction Machine StanVac 1 at wholesale competitive prices from India.

### **Specifications :**

Portable Slow Suction Machine StanVac 1 is made of mild steel with electrostatically powder coated mounted on Rubber pad

Portable Slow Suction Machine StanVac 1 has Oil Immersed Motorised Noiseless Pump

Portable Slow Suction Machine StanVac 1 is fitted with Reuseable Filter Jars : 2 Jars.

The Polycarbonate Jars are fitted on the top of the Mild Steel Base with lids with over flow safety device.

Vacuum Gauge has 50 mm diameter graduated in mm Hg.

Vacuum is -700 mm Hg  $\pm$  10% which is controlled by knob.

The Suction Capacity is 20 Ltr/min

Power : 230V-50 Hz

We are Manufacturer, Exporters and Suppliers of Portable Slow Suction Machine StanVac 1 at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)