

## Labycare-Portable Slow Suction Machine DeVac (SKSS-SM27)



Portable Slow Suction Machine DeVac Manufacturer by LABYCARE

### **Description :**

Portable Slow Suction Machine DeVac is made by Fibre Reinforced Plastic (FRP) shock proof mounted on heavy duty 65mm castors

Portable Slow Suction Machine DeVac is fitted with Reuseable Filter Jars : 2 x 2 Ltrs.

The Polycarbonate Jar is fitted with lids on the base of the Suction Unit with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 100 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Rocker Piston Pump

This machine is ideal for MTP / medical / surgical procedures

Vacuum Capacity is -710 mm Hg  $\pm$  10% at 30-35 LPM which is controlled by knob.

Noise Level < 50 dB A  $\pm$  3 Almost Whispers

Power : 220/230V , 50Hz, 1 Phase

We are Manufacturer, Exporters and Suppliers of Portable Slow Suction Machine DeVac at wholesale competitive prices from India.

### **Specifications :**

Portable Slow Suction Machine DeVac is made by Fibre Reinforced Plastic (FRP) shock proof mounted on heavy duty 65mm castors

Portable Slow Suction Machine DeVac is fitted with Reuseable Filter Jars : 2 x 2 Ltrs.

The Polycarbonate Jar is fitted with lids on the base of the Suction Unit with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 100 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Rocker Piston Pump

This machine is ideal for MTP / medical / surgical procedures

Vacuum Capacity is -710 mm Hg  $\pm$  10% at 30-35 LPM which is controlled by knob.

Noise Level < 50 dB A  $\pm$  3 Almost Whispers

Power : 220/230V , 50Hz, 1 Phase

We are Manufacturer, Exporters and Suppliers of Portable Slow Suction Machine DeVac at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)