

## Labycare-Foot Suction MS Base Glass Jar (SKSS-SM28)



Foot Suction MS Base Glass Jar Manufacturer by LABYCARE

### **Description :**

Foot Suction Machine is manually operated

Foot Suction Base is Mild Steel Powder Coated

Suction Capacity is -600mmHg  $\pm$  10 at 200ml per stroke

Foot Suction is fitted with Piston Pump

Foot Suction comes with Glass Jars 2 x 600ml capacity with over flow system & ABS lid

Tubing is non collapsible

Vacuum Gauge has 50 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

We are Manufacturer, Exporters and Suppliers of Foot Suction MS Base Glass Jar at wholesale competitive prices from India.

### **Specifications :**

Foot Suction Machine is manually operated

Foot Suction Base is Mild Steel Powder Coated

Suction Capacity is -600mmHg  $\pm$  10 at 200ml per stroke

Foot Suction is fitted with Piston Pump

Foot Suction comes with Glass Jars 2 x 600ml capacity with over flow system & ABS lid

Tubing is non collapsible

Vacuum Gauge has 50 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

We are Manufacturer, Exporters and Suppliers of Foot Suction MS Base Glass Jar at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)