

Labycare-Foot Suction MS Base PC Jar (SKSS-SM29)



Foot Suction MS Base PC Jar Manufacturer by LABYCARE

Description :

Foot Suction Machine is manually operated

Foot Suction Base is Mild Steel Powder Coated

Suction Capacity is -600mmHg \pm 10 at 200ml per stroke

Foot Suction is fitted with Piston Pump

Foot Suction comes with Polycarbonate Jars 2 x 600ml capacity with over flow system & ABS lid

Tubing is non collapsible

Vacuum Gauge has 50 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

We are Manufacturer, Exporters and Suppliers of Foot Suction MS Base PC Jar at wholesale competitive prices from India.

Specifications :

Foot Suction Machine is manually operated

Foot Suction Base is Mild Steel Powder Coated

Suction Capacity is -600mmHg \pm 10 at 200ml per stroke

Foot Suction is fitted with Piston Pump

Foot Suction comes with Polycarbonate Jars 2 x 600ml capacity with over flow system & ABS lid

Tubing is non collapsible

Vacuum Gauge has 50 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

We are Manufacturer, Exporters and Suppliers of Foot Suction MS Base PC Jar at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com