

## Labycare-Foot Suction EcoVac (SKSS-SM32)



Foot Suction EcoVac Manufacturer by LABYCARE

### **Description :**

Foot Suction EcoVac is manually operated

Foot Suction EcoVac Base is Heavy Mild Steel Powder Coated

Suction Capacity is -600mmHg  $\pm$  10 at 200ml per stroke

Foot Suction EcoVac is fitted with heavy Piston Pump

Foot Suction EcoVac comes with Polycarbonate Jar 1 x 1500ml capacity with over flow system & ABS lid

Tubing is non collapsible

Vacuum Gauge has 50 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

We are Manufacturer, Exporters and Suppliers of Foot Suction EcoVac at wholesale competitive prices from India.

### **Specifications :**

Foot Suction EcoVac is manually operated

Foot Suction EcoVac Base is Heavy Mild Steel Powder Coated

Suction Capacity is -600mmHg  $\pm$  10 at 200ml per stroke

Foot Suction EcoVac is fitted with heavy Piston Pump

Foot Suction EcoVac comes with Polycarbonate Jar 1 x 1500ml capacity with over flow system & ABS lid

Tubing is non collapsible

Vacuum Gauge has 50 mm diameter graduated from 0 - 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

We are Manufacturer, Exporters and Suppliers of Foot Suction EcoVac at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)