

Labycare-Ward Vacuum Unit (SKSS-SM34)



Ward Vacuum Unit Manufacturer by LABYCARE

Description :

Ward Vacuum Unit is made of Polycarbonate

Ward Vacuum Unit is fitted with Regulator which is made of Aluminium & ABS Plastic

Vacuum is controlled by Regulator Knob

Vacuum Gauge has 50 mm diameter graduated in mmHg.

Ward Vacuum Unit is used in Pipe Line Fitting

We are Manufacturer, Exporters and Suppliers of Ward Vacuum Unit at wholesale competitive prices from India.

Specifications :

Ward Vacuum Unit is made of Polycarbonate

Ward Vacuum Unit is fitted with Regulator which is made of Aluminium & ABS Plastic

Vacuum is controlled by Regulator Knob

Vacuum Gauge has 50 mm diameter graduated in mmHg.

Ward Vacuum Unit is used in Pipe Line Fitting

We are Manufacturer, Exporters and Suppliers of Ward Vacuum Unit at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com