

Labycare-Vacuum Extractor Manual SS Cups (SKSS-SM35)



Vacuum Extractor Manual SS Cups Manufacturer by LABYCARE

Description :

Vacuum Extractor Manual SS Cups is based on manual operation.

Vacuum Extractor is supplied with three Stainless Steel Cups

Glass Jar is fitted with vacuum gauge and control valve

Stand for Glass Jar is provided.

The Size of Cups Are 40mm, 50mm & 60mm

Vacuum Pump is Manually Operated

Vacuum Pump and Bottles is fitted with Rubber tubes

Vacuum Extractor has autoclavable glass jar with Vacuum Gauge

Vacuum Gauge is Bourden type has 50 mm diameter graduated from 0 - 760 mmHg 0-100kPa

The complete unit is packed in housing case

We are Manufacturer, Exporters and Suppliers of Vacuum Extractor Manual SS Cups at wholesale competitive prices from India.

Specifications :

Vacuum Extractor Manual SS Cups is based on manual operation.

Vacuum Extractor is supplied with three Stainless Steel Cups

Glass Jar is fitted with vacuum gauge and control valve

Stand for Glass Jar is provided.

The Size of Cups Are 40mm, 50mm & 60mm

Vacuum Pump is Manually Operated

Vacuum Pump and Bottles is fitted with Rubber tubes

Vacuum Extractor has autoclavable glass jar with Vacuum Gauge

Vacuum Gauge is Bourden type has 50 mm diameter graduated from 0 - 760 mmHg 0-100kPa

The complete unit is packed in housing case

We are Manufacturer, Exporters and Suppliers of Vacuum Extractor Manual SS Cups at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com