

Labycare-Vacuum Extractor Manual Silicone Cups (SKSS-SM36)



Vacuum Extractor Manual Silicone Cups Manufacturer by LABYCARE

Description :

Vacuum Extractor Manual is based on manual operation.

Vacuum Extractor is supplied with three Silicone Cups.

Glass Jar is fitted with vacuum gauge and control valve

Stand for Glass Jar is provided.

The Size of Cups Are 40mm, 50mm & 60mm.

Vacuum Pump is Manually Operated.

Vacuum Pump and Bottles is fitted with Rubber tubes.

Vacuum Extractor has autoclavable glass jar with Vacuum Gauge.

Vacuum Gauge is Bourden type has 50 mm diameter graduated from 0 - 760 mmHg 0-100kPa.

The complete unit is packed in housing case.

We are Manufacturer, Exporters and Suppliers of Vacuum Extractor Manual Silicone Cups at wholesale competitive prices from India.

Specifications :

Vacuum Extractor Manual is based on manual operation.

Vacuum Extractor is supplied with three Silicone Cups.

Glass Jar is fitted with vacuum gauge and control valve

Stand for Glass Jar is provided.

The Size of Cups Are 40mm, 50mm & 60mm.

Vacuum Pump is Manually Operated.

Vacuum Pump and Bottles is fitted with Rubber tubes.

Vacuum Extractor has autoclavable glass jar with Vacuum Gauge.

Vacuum Gauge is Bourden type has 50 mm diameter graduated from 0 - 760 mmHg 0-100kPa.

The complete unit is packed in housing case.

We are Manufacturer, Exporters and Suppliers of Vacuum Extractor Manual Silicone Cups at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com