

## Labycare-Vacuum Extractor Manual Bird Type SS Cups (SKSS-SM37)



Vacuum Extractor Manual Bird Type SS Cups Manufacturer by LABYCARE

### **Description :**

Vacuum Extractor Set is based on manual operation.

Vacuum Extractor is supplied with three Stainless Steel Cups

The Size of Cups Are 50mm, 60mm & 70mm

Vacuum Pump is Manually Operated

Vacuum Extractor has autoclavable glass jar with Vacuum Gauge

Vacuum Gauge is Bourden type has 50 mm diameter graduated from 0 - 760 mmHg 0-100kPa

All parts coming into contact with secretion are autoclave safe.

The complete unit is packed in Superior Quality housing case

We are Manufacturer, Exporters and Suppliers of Vacuum Extractor Manual Bird Type SS Cups at wholesale competitive prices from India.

### **Specifications :**

Vacuum Extractor Set is based on manual operation.

Vacuum Extractor is supplied with three Stainless Steel Cups

The Size of Cups Are 50mm, 60mm & 70mm

Vacuum Pump is Manually Operated

Vacuum Extractor has autoclavable glass jar with Vacuum Gauge

Vacuum Gauge is Bourden type has 50 mm diameter graduated from 0 - 760 mmHg 0-100kPa

All parts coming into contact with secretion are autoclave safe.

The complet unit is packed in Superior Quality housing case

We are Manufacturer, Exporters and Suppliers of Vacuum Extractor Manual Bird Type SS Cups at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)