

## Labycare-Vertical Foot Suction Pump StanVac (SKSS-SM46)



Vertical Foot Suction Pump StanVac Manufacturer by LABYCARE

### **Description :**

Vertical Foot Suction Pump StanVac is made by mild steel with electrostatically powder coated

Vertical Foot Suction Pump StanVac is fitted with Reuseable Filter Jars : 1 x 1.6 Ltrs.

The Jar is fitted with synthetic rubber lids with over flow safety device.

Vacuum Gauge has 50 mm diameter graduated in mm Hg.

Vertical Foot Suction Pump comes with long pressure tubing with synthetic rubber lids.

We are Manufacturer, Exporters and Suppliers of Vertical Foot Suction Pump StanVac at wholesale competitive prices from India.

### **Specifications :**

Vertical Foot Suction Pump StanVac is made by mild steel with electrostatically powder coated

Vertical Foot Suction Pump StanVac is fitted with Reuseable Filter Jars : 1 x 1.6 Ltrs.

The Jar is fitted with synthetic rubber lids with over flow safety device.

Vacuum Gauge has 50 mm diameter graduated in mm Hg.

Vertical Foot Suction Pump comes with long pressure tubing with synthetic rubber lids.

We are Manufacturer, Exporters and Suppliers of Vertical Foot Suction Pump StanVac at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)