

Labycare-Vertical Foot Suction Pump DeVac (SKSS-SM47)



Vertical Foot Suction Pump DeVac Manufacturer by LABYCARE

Description :

Vertical Foot Suction Pump DeVac Made of Stainless Steel base

Vertical Foot Suction Pump DeVac is fitted with Reuseable Filter Jars : 1 x 1.6 Ltrs.

The Jar is fitted with synthetic rubber lids with over flow safety device.

Vacuum Gauge has 50 mm diameter graduated in mm Hg.

Vertical Foot Suction Pump comes with long pressure tubing with lids.

We are Manufacturer, Exporters and Suppliers of Vertical Foot Suction Pump DeVac at wholesale competitive prices from India.

Specifications :

Vertical Foot Suction Pump DeVac Made of Stainless Steel base

Vertical Foot Suction Pump DeVac is fitted with Reuseable Filter Jars : 1 x 1.6 Ltrs.

The Jar is fitted with synthetic rubber lids with over flow safety device.

Vacuum Gauge has 50 mm diameter graduated in mm Hg.

Vertical Foot Suction Pump comes with long pressure tubing with lids.

We are Manufacturer, Exporters and Suppliers of Vertical Foot Suction Pump DeVac at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com