

Labycare-Hand Held Suction Device StanVac (SKSS-SM48)



Hand Held Suction Device StanVac Manufacturer by LABYCARE

Description :

Hand Held Suction Device is made by superior quality Plastic.

Hand Held Suction Device StanVac is sturdy and Non tiring unit.

Hand held Suction Unit is designed to provide immediate, quick and effective solution in emergency.

Volume of Disposable Container is 300 ML.

Operating Environmental Temperature is -20 degrees C to +50 degrees C

Storage Environmental temperature is -40o C to +60o C

50% For Child, 100% For Adult,

Child Catheter connection External Diameter ø 9 mm Length 21cm

Adult Catheter connection External Diameter ø 17 mm Length 21cm

We are Manufacturer, Exporters and Suppliers of Hand Held Suction Device StanVac at wholesale competitive prices from India.

Specifications :

Hand Held Suction Device is made by superior quality Plastic.

Hand Held Suction Device StanVac is sturdy and Non tiring unit.

Hand held Suction Unit is designed to provide immediate, quick and effective solution in emergency.

Volume of Disposable Container is 300 ML.

Operating Environmental Temperature is -20 degrees C to +50 degrees C

Storage Environmental temperature is -40o C to +60o C

50% For Child, 100% For Adult,

Child Catheter connection External Diameter \varnothing 9 mm Length 21cm

Adult Catheter connection External Diameter \varnothing 17 mm Length 21cm

We are Manufacturer, Exporters and Suppliers of Hand Held Suction Device StanVac at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com