

Labycare-Suction Foot Switch (SKSS-SM55)



Suction Foot Switch

Description :

Suction Foot Switch is made of Shockproof ABS Plastic

Suction Foot Switch is fitted with 10A Micro switch

The Length of wire is 2 meters

Suction Foot Switch is mounted on rubber base.

We are Manufacturer, Exporters and Suppliers of Suction Foot Switch at wholesale competitive prices from India.

Specifications :

Suction Foot Switch is made of Shockproof ABS Plastic

Suction Foot Switch is fitted with 10A Micro switch

The Length of wire is 2 meters

Suction Foot Switch is mounted on rubber base.

We are Manufacturer, Exporters and Suppliers of Suction Foot Switch at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com