

Labycare-Autopsy Tables Electric (SKSSATE378)



Autopsy table electric with down-draft and height adjustment. Table completely made of stainless steel 1.4301 (AISI 304). Optionally it can be made of stainless steel 1.4404 (AISI 316). The table features a roller system where the body tray can be directly loaded without need to lift the body manually, specially designed to increase the ergonomics when transferring bodies from the transport trolley to the autopsy table. Under the roller system, the table has a slope, one outflow and integrated sprinkler system. At the end of the table, basin with cold water tap and overflow pipe. Under the basin, there is a mixing tap for cold and hot water with thermostat and a hand sprinker with 3 m shower plastic hose. The height adjustment is reached by an electric motor operated by manual buttons on the side of the column. On the central support column 2 waterproof electric sockets. It is possible to connect the down draft system to the building ventilation system.

Description :

utopsy table electric with down-draft and height adjustment. Table completely made of stainless steel 1.4301). Optionally it can be made of stainless steel 1.4404). The table features a roller system where the body tray can be directly loaded without need to lift the body manually, specially designed to increase the ergonomics when transferring bodies from the transport trolley to the autopsy table. Under the roller system, the table has a slope, one outflow and integrated sprinkler system. At the end of the table, basin with cold water tap and overflow pipe. Under the basin, there is a mixing tap for cold and hot water with thermostat and a hand sprinker with 3 m shower plastic hose. The height adjustment is reached by an electric motor operated by manual buttons on the side of the column. On the central support column 2 waterproof electric sockets. It is possible to connect the down draft system to the building ventilation system.

Specifications :

--	--

Portable	Yes
Material	Stainless
Brand	Weswox
Table Length	91 Inch
Width	30 Inch
Height	36 Inch
Material Grade	SS 304

DISCRIPTION:

-All S.S. Construction sturdy support framework. One endmounted intergrated large SS wash basin with one free rotating faucet and one long reach hand shower. Built in under table plumbing and drainage. (External pluming not included).Comes with one SS tubular portable and adjustabel overhead light on casters.