

Labycare-Analytical Balances with Internal Calibration (SKSS-HS01)



Analytical Balances with Internal Calibration Manufacturer by LABYCARE

Description :

Analytical Balance is also known as Laboratory Digital Balance

Automatic internal calibration with built in weight

It is LCD Display with back light

The Standard RS 232 C Interface is available in this model

This is based on E.M.F.C. Technology

The balance has S.S weighing pan, Level Indicator, adjustable fit

Full tare range up to maximum capacity of the Balance

Repeatability : $\pm 0.1\text{mg}$

Linearity : $\pm 0.2\text{ mg}$

Pan size : 85ømm

Stabilization Time : 3-4 Sec

Communication : RS 232 Interface

Calibration : Internal calibration

Weight Units : Multi unit

Power : 230 V / 11V AC or 120V/ 11V AC

Display : LCD Display with back light

Working Temperature : 10° C - 40° C

We are Manufacturer, Exporters and Suppliers of Analytical Balances with Internal Calibration at wholesale competitive prices from India.

Specifications :

Analytical Balance is also known as Laboratory Digital Balance

Automatic internal calibration with built in weight

It is LCD Display with back light

The Standard RS 232 C Interface is available in this model

This is based on E.M.F.C. Technology

The balance has S.S weighing pan, Level Indicator, adjustable fit

Full tare range up to maximum capacity of the Balance

Repeatability : $\pm 0.1\text{mg}$

Linearity : $\pm 0.2\text{ mg}$

Pan size : 85ømm

Stabilization Time : 3-4 Sec

Communication : RS 232 Interface

Calibration : Internal calibration

Weight Units : Multi unit

Power : 230 V / 11V AC or 120V/ 11V AC

Display : LCD Display with back light

Working Temperature : 10° C - 40° C

We are Manufacturer, Exporters and Suppliers of Analytical Balances with Internal Calibration at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com