

Labycare-Baby Scale Digital 20 Kg 10 Gms (SKSS-HS09)



Baby Scale Digital 20 Kg 10 Gms

Description :

Baby Scale Digital is made of ABS Plastic

Baby Scale Digital is also known as Digital Baby Weighing Scale

It is Electronic Digital Non Slip which can weight child upto 20 Kg

This Weighing Scale is equipped with a high precision strain gauge and has high performance CPU ensuring accurate weight.

The scale features on/off button.

It has the ability to display weight in kilograms, pounds or stones

It comes with an overload indicator and a low power indicator.

It is easy to clean.

Capacity: 20Kg/44lbs

Auto off function

Low-battery indication

Overload protective system

Weighing Units: Kg/Lbs/St

Division: 10g

Battery: 2 x AAA

We are Manufacturer, Exporters and Suppliers of Baby Scale Digital 20 Kg 10 Gms at wholesale competitive prices from India.

Specifications :

Baby Scale Digital is made of ABS Plastic

Baby Scale Digital is also known as Digital Baby Weighing Scale

It is Electronic Digital Non Slip which can weight child upto 20 Kg

This Weighing Scale is equipped with a high precision strain gauge and has high performance CPU ensuring accurate weight.

The scale features on/off button.

It has the ability to display weight in kilograms, pounds or stones

It comes with an overload indicator and a low power indicator.

It is easy to clean.

Capacity: 20Kg/44lbs

Auto off function

Low-battery indication

Overload protective system

Weighing Units: Kg/Lbs/St

Division: 10g

Battery: 2 x AAA

We are Manufacturer, Exporters and Suppliers of Baby Scale Digital 20 Kg 10 Gms at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com