

Labycare-Height Measuring Scale Digital BMI (SKSS-HS43)



Height Measuring Scale Digital BMI Manufacturer by LABYCARE

Description :

Adult digital scale with height and weight

- Equipped with high precision strain gauge sensor system
- Built-in 8mm tempered safety glass
- Capacity: 3kg-200kg (scale in the state, it can be 0.3kg-200kg)
- Division of weight: 50g
- Height measuring range: 0.80-1.95 meter (division of height: 5mm)
- LCD display Size: 105*69.5 mm (back light)
- Height measuring function
- Body fat measuring and analysing function
- Voice function
- Single set of data recording/checking function
- Automatic on-step, continuous weighing Auto-zero
- Low battery/over-load indicator
- Power Supply- Power by 4*AAA alkaline batteries or 6V•15V DC power supply. Before using, please open the battery cover to fix battery or plug in DC power."

We are Manufacturer, Exporters and Suppliers of Height Measuring Scale Digital BMI at wholesale competitive prices from India.

Specifications :

Adult digital scale with height and weight

- Equipped with high precision strain gauge sensor system
- Built-in 8mm tempered safety glass
- Capacity: 3kg-200kg (scale in the state, it can be 0.3kg-200kg)
- Division of weight: 50g
- Height measuring range: 0.80-1.95 meter (division of height: 5mm)
- LCD display Size: 105*69.5 mm (back light)
- Height measuring function
- Body fat measuring and analysing function
- Voice function
- Single set of data recording/checking function
- Automatic on-step, continuous weighing Auto-zero
- Low battery/over-load indicator
- Power Supply- Power by 4*AAA alkaline batteries or 6V•15V DC power supply. Before using, please open the battery cover to fix battery or plug in DC power."

We are Manufacturer, Exporters and Suppliers of Height Measuring Scale Digital BMI at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com