

Labycare-Weighing Scales Mechanical Column Type (SKSS-HS44)



Weighing Scales Mechanical Column Type Manufacturer by LABYCARE

Description :

Weighing Scales Mechanical Column Type two in one weighing and height scale

The scale has the advantages of large weighing range, accuracy and high sensitivity.

The health scale is also equipped with pulley which can make the movement in short distance more convenient.

Max.weighing:160-200kg

Scale division(d) & Accuracy :100g

Height range to be measured:75-200cm

Min.value of height per division:0.5cm

Size of packing box:1500×510×155mm

Load platform area(L×W):375×275mm

Measurement(L×W×H):530×275×1485mm

We are Manufacturer, Exporters and Suppliers of Weighing Scales Mechanical Column Type at wholesale competitive prices from India.

Specifications :

Weighing Scales Mechanical Column Type two in one weighing and height scale

The scale has the advantages of large weighing range, accuracy and high sensitivity.

The health scale is also equipped with pulley which can make the movement in short distance more convenient.

Max. weighing: 160-200kg

Scale division(d) & Accuracy : 100g

Height range to be measured: 75-200cm

Min. value of height per division: 0.5cm

Size of packing box: 1500×510×155mm

Load platform area(L×W): 375×275mm

Measurement(L×W×H): 530×275×1485mm

We are Manufacturer, Exporters and Suppliers of Weighing Scales Mechanical Column Type at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com