

Labycare-Weighing Scales Mechanical Dial Type (SKSS-HS45)



Weighing Scales Mechanical Dial Type Manufacturer by LABYCARE

Description :

Weighing Scales Mechanical Dial Type two in one weighing and height scale

The scale has the advantages of large weighing range, accuracy and high sensitivity.

Max.weighing:120/150/160kg

Min.value per division:0.5kg

Aluminous Height range to be measured:80-210cm

Min. Value of height per division:0.5cm

Load platform area(L×W):375×275mm

Measurement(L×W×H):695×286×930mm

We are Manufacturer, Exporters and Suppliers of Weighing Scales Mechanical Dial Type at wholesale competitive prices from India.

Specifications :

Weighing Scales Mechanical Dial Type two in one weighing and height scale

The scale has the advantages of large weighing range, accuracy and high sensitivity.

Max.weighing:120/150/160kg

Min.value per division:0.5kg

Aluminous Height range to be measured:80-210cm

Min. Value of height per division:0.5cm

Load platform area(L×W):375×275mm

Measurement(L×W×H):695×286×930mm

We are Manufacturer, Exporters and Suppliers of Weighing Scales Mechanical Dial Type at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com